


Tos-PEG2-OH

Chemical Properties

CAS No. :	118591-58-5
Formula:	C11H16O5S
Molecular Weight:	260.31
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Tos-PEG2-OH, a PEG-based linker for PROTACs, connects two essential ligands critical for forming PROTAC molecules, enabling selective protein degradation through the ubiquitin-proteasome system [within cells].
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two distinct ligands connected by a linker: one ligand targets an E3 ubiquitin ligase, while the other binds to the target protein. By leveraging the intracellular ubiquitin-proteasome system, PROTACs facilitate the selective degradation of target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8416 mL	19.2079 mL	38.4157 mL
5 mM	0.7683 mL	3.8416 mL	7.6831 mL
10 mM	0.3842 mL	1.9208 mL	3.8416 mL
50 mM	0.0768 mL	0.3842 mL	0.7683 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481