


Vinagensenoside R8

Chemical Properties

CAS No. : 156042-22-7
 Formula: C₄₈H₈₂O₁₉
 Molecular Weight: 963.165
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Vinagensenoside R8 is a natural product
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0382 mL	5.1912 mL	10.3824 mL
5 mM	0.2076 mL	1.0382 mL	2.0765 mL
10 mM	0.1038 mL	0.5191 mL	1.0382 mL
50 mM	0.0208 mL	0.1038 mL	0.2076 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481