


SOS1-IN-7

Chemical Properties

CAS No. : 2755415-30-4
 Formula: C₂₃H₂₅F₃N₄O₃
 Molecular Weight: 462.46
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	SOS1-IN-7 (compound 18-p1) is a potent SOS1 inhibitor with IC ₅₀ values of 20 nM for SOS1-G12D and 67 nM for SOS1-G12V.
Targets(IC ₅₀)	Others,Ras

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1623 mL	10.8117 mL	21.6235 mL
5 mM	0.4325 mL	2.1623 mL	4.3247 mL
10 mM	0.2162 mL	1.0812 mL	2.1623 mL
50 mM	0.0432 mL	0.2162 mL	0.4325 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481