

CCG-273220

## Chemical Properties


CAS No. : 2750414-35-6

Formula: C<sub>26</sub>H<sub>24</sub>ClFN<sub>4</sub>O<sub>3</sub>

Molecular Weight: 494.95

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	CCG-273220 is a covalent inhibitor of G protein-coupled receptor (GPCR) kinase 5 (GRK5) (IC <sub>50</sub> : 220 nM). CCG-273220 covalently binds the GRK5 subfamily-specific residue Cys474 and is more selective for GRK5 than GRK2.
Targets(IC <sub>50</sub> )	Others,GRK

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0204 mL	10.102 mL	20.2041 mL
5 mM	0.4041 mL	2.0204 mL	4.0408 mL
10 mM	0.202 mL	1.0102 mL	2.0204 mL
50 mM	0.0404 mL	0.202 mL	0.4041 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481