


Tert-Pentylbenzene

Chemical Properties

CAS No. :	2049-95-8
Formula:	C ₁₁ H ₁₆
Molecular Weight:	148.25
Storage:	Pure form: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Tert-Pentylbenzene is a benzene derivative widely used in biochemical experiments and drug synthesis research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.7454 mL	33.7268 mL	67.4536 mL
5 mM	1.3491 mL	6.7454 mL	13.4907 mL
10 mM	0.6745 mL	3.3727 mL	6.7454 mL
50 mM	0.1349 mL	0.6745 mL	1.3491 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481