

cis-3-Hexenyl acetate

Chemical Properties


CAS No. : 3681-71-8

Formula: C₈H₁₄O₂

Molecular Weight: 142.20

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	cis-3-Hexenyl acetate is an unsaturated ester compound widely used in biochemical experiments and drug synthesis research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.0323 mL	35.1617 mL	70.3235 mL
5 mM	1.4065 mL	7.0323 mL	14.0647 mL
10 mM	0.7032 mL	3.5162 mL	7.0323 mL
50 mM	0.1406 mL	0.7032 mL	1.4065 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Mairean CP, et al. Gas-Phase Kinetics of a Series of cis-3-Hexenyl Esters with OH Radicals under Simulated Atmospheric Conditions. J Phys Chem A. 2024 Aug 1;128(30):6274-6285.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481