

2-Isopropyl-6-methylaniline

Chemical Properties

CAS No. :	5266-85-3
Formula:	C10H15N
Molecular Weight:	149.233
Storage:	Pure form: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.

Biological Description

Description	2-Isopropyl-6-methylaniline is an aromatic amine compound widely used in biochemical experiments and drug synthesis research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.7011 mL	33.5053 mL	67.0107 mL
5 mM	1.3402 mL	6.7011 mL	13.4021 mL
10 mM	0.6701 mL	3.3505 mL	6.7011 mL
50 mM	0.134 mL	0.6701 mL	1.3402 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481