


2-Phenoxyphenol

Chemical Properties

CAS No. :	2417-10-9
Formula:	C ₁₂ H ₁₀ O ₂
Molecular Weight:	186.21
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Phenoxyphenol has antibacterial activity and inhibitory activity against Staphylococcus aureus, Staphylococcus aureus, Escherichia coli, Burkholderia pseudomallei, etc. It is widely used in biochemical experiments and drug synthesis research.
Targets(IC50)	Antibacterial

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.3703 mL	26.8514 mL	53.7028 mL
5 mM	1.0741 mL	5.3703 mL	10.7406 mL
10 mM	0.537 mL	2.6851 mL	5.3703 mL
50 mM	0.1074 mL	0.537 mL	1.0741 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Yu B, et al. Selected dechlorination of triclosan by high-performance g-C₃N₄/Bi₂MoO₆ composites: Mechanisms and pathways. Chemosphere. 2023 Jan;312(Pt 1):137247.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481