


## Benzo[c]phenanthrene

## Chemical Properties

CAS No. :	195-19-7
Formula:	C <sub>18</sub> H <sub>12</sub>
Molecular Weight:	228.29
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Benzo[c]phenanthrene exhibits carcinogenic activity through skin contact and is widely used in biochemical experiments and drug synthesis research.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.3804 mL	21.902 mL	43.8039 mL
5 mM	0.8761 mL	4.3804 mL	8.7608 mL
10 mM	0.438 mL	2.1902 mL	4.3804 mL
50 mM	0.0876 mL	0.438 mL	0.8761 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

FLURRY R Jr. A SIMPLE METHOD FOR PREDICTING THE CARCINOGENIC PROPERTIES OF POLYCYCLIC AROMATIC MOLECULES. J Med Chem. 1964 Sep;7:668-70.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481