


2-Methyl-1-indanone

Chemical Properties

CAS No. :	17496-14-9
Formula:	C ₁₀ H ₁₀ O
Molecular Weight:	146.19
Storage:	Store at 4°C <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	2-Methyl-1-indanone is an acetylcholinesterase inhibitor and can be used for biochemical experiments and drug synthesis.
Targets(IC50)	Cholinesterase (ChE)

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.8404 mL	34.2021 mL	68.4041 mL
5 mM	1.3681 mL	6.8404 mL	13.6808 mL
10 mM	0.684 mL	3.4202 mL	6.8404 mL
50 mM	0.1368 mL	0.684 mL	1.3681 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Cardozo MG, et al. Quantitative structure-activity relationship. QSAR, analyses of the substituted indanone and benzylpiperidine rings of a series of indanone-benzylpiperidine inhibitors of acetylcholinesterase. Journal of Medicinal Chemistry. 1992 Feb;35(3):584-589.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481