

5H-Pyrido[4,3-b]indole

Chemical Properties

CAS No. :	244-69-9
Formula:	C ₁₁ H ₈ N ₂
Molecular Weight:	168.195
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	5H-Pyrido[4,3-b]indole is a potential AchE/ChE inhibitor with potential antiviral activity for the study of neurodegenerative diseases.
Targets(IC50)	Cholinesterase (ChE)

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.9453 mL	29.7265 mL	59.453 mL
5 mM	1.1891 mL	5.9453 mL	11.8906 mL
10 mM	0.5945 mL	2.9727 mL	5.9453 mL
50 mM	0.1189 mL	0.5945 mL	1.1891 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481