

2,5-Diphenylpyridine

Chemical Properties

CAS No. :	15827-72-2
Formula:	C17H13N
Molecular Weight:	231.292
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.

Biological Description

Description	2,5-Diphenylpyridine has anticancer activity against a variety of tumour cell lines when used in combination with platinum, and is widely used in biochemical experiments and drug synthesis research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.3236 mL	21.6179 mL	43.2358 mL
5 mM	0.8647 mL	4.3236 mL	8.6472 mL
10 mM	0.4324 mL	2.1618 mL	4.3236 mL
50 mM	0.0865 mL	0.4324 mL	0.8647 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Zou G, et al. Columnar Iridium(III) Metallomesogens Based on Polycatenar Pyridyltetrazolate with Ambipolar Carrier Mobility Behavior. *Inorg Chem.* 2019 Jan 7;58(1):861-869.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481