

## 4-(Pyridin-2-yl)butan-1-ol

## Chemical Properties


CAS No. : 17945-79-8

Formula: C<sub>9</sub>H<sub>13</sub>NO

Molecular Weight: 151.206

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	4-(Pyridin-2-yl)butan-1-ol is a pyridine derivative widely used in biochemical experiments and drug synthesis research.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.6133 mL	33.0666 mL	66.1332 mL
5 mM	1.3227 mL	6.6133 mL	13.2266 mL
10 mM	0.6613 mL	3.3067 mL	6.6133 mL
50 mM	0.1323 mL	0.6613 mL	1.3227 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481