


## 9,10-Dimethylantracene

## Chemical Properties

CAS No. :	781-43-1
Formula:	C <sub>16</sub> H <sub>14</sub>
Molecular Weight:	206.28
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	9,10-Dimethylantracene has antibacterial activity against Escherichia coli and is widely used in biochemical experiments and drug synthesis research.
Targets(IC50)	Antibacterial

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8478 mL	24.2389 mL	48.4778 mL
5 mM	0.9696 mL	4.8478 mL	9.6956 mL
10 mM	0.4848 mL	2.4239 mL	4.8478 mL
50 mM	0.097 mL	0.4848 mL	0.9696 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Albiter E, et al. Photosensitized oxidation of 9,10-dimethylantracene with singlet oxygen by using a safranin O/silica composite under visible light. Photochem Photobiol Sci. 2015 Mar; 14(3):597-602.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481