

## 5-Ethylindole

## Chemical Properties


CAS No. : 68742-28-9

Formula: C<sub>10</sub>H<sub>11</sub>N

Molecular Weight: 145.2

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	5-Ethylindole is an indole derivative widely used in biochemical experiments and drug synthesis research.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.8871 mL	34.4353 mL	68.8705 mL
5 mM	1.3774 mL	6.8871 mL	13.7741 mL
10 mM	0.6887 mL	3.4435 mL	6.8871 mL
50 mM	0.1377 mL	0.6887 mL	1.3774 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481