

## Anti-VZV-gE Antibody (93K)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:   
Store at low temperature  
-20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-VZV-gE Antibody (93K) is an antibody. Its expression system is HEK293 cells, the clone number is 93K, and the isotype control is Mouse & Human IgG1 Chimera.
Targets(IC50)	HSV

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481