

## Anti-HMPV-F Antibody (MPV458)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-HMPV-F Antibody (MPV458) is an antibody. Its expression system is CHO cells, the clone number is MPV458, and the isotype control is IgG1, Kappa.
Targets(IC50)	Antiviral

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481