

## Anti-HMPV-F Antibody (MPV467)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-HMPV-F Antibody (MPV467) is an antibody that specifically targets the human HMPV-F protein and is capable of inhibiting HMPV-F replication in the lungs in animal models.
Targets(IC50)	Antiviral

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481