

Anti-HMPV-F Antibody (MPV487)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-HMPV-F Antibody (MPV487) is an antibody. Its expression system is CHO cells, the clone number is MPV487, and the isotype control is IgG1, Kappa.
Targets(IC50)	Antiviral

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481