

## Cinnamyl 3-methylbutanoate

## Chemical Properties


CAS No. : 140-27-2

Formula: C<sub>14</sub>H<sub>18</sub>O<sub>2</sub>

Molecular Weight: 218.29

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Cinnamyl 3-methylbutanoate is an unsaturated ester widely used in biochemical experiments and drug synthesis research.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.5811 mL	22.9053 mL	45.8106 mL
5 mM	0.9162 mL	4.5811 mL	9.1621 mL
10 mM	0.4581 mL	2.2905 mL	4.5811 mL
50 mM	0.0916 mL	0.4581 mL	0.9162 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481