


## D-(+)-Cellotetraose tetradecaacetate

### Chemical Properties

CAS No. : 83058-25-7  
 Formula: C<sub>52</sub>H<sub>70</sub>O<sub>35</sub>  
 Molecular Weight: 1255.09  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



### Biological Description

Description	D-(+)-Cellotetraose tetradecaacetate is a reagent used in biochemical reactions.
Targets(IC50)	Others

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7968 mL	3.9838 mL	7.9676 mL
5 mM	0.1594 mL	0.7968 mL	1.5935 mL
10 mM	0.0797 mL	0.3984 mL	0.7968 mL
50 mM	0.0159 mL	0.0797 mL	0.1594 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481