


L-Idaric acid, 1,4-lactone

Chemical Properties

CAS No. :	80876-59-1
Formula:	C ₆ H ₈ O ₇
Molecular Weight:	192.12
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	L-Idaric acid,1,4-lactone is a reagent used in biochemical reactions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2051 mL	26.0254 mL	52.0508 mL
5 mM	1.041 mL	5.2051 mL	10.4102 mL
10 mM	0.5205 mL	2.6025 mL	5.2051 mL
50 mM	0.1041 mL	0.5205 mL	1.041 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481