

5-Keto-D-gluconic acid potassium

Chemical Properties


CAS No. : 91446-96-7

Formula: C₆H₉KO₇

Molecular Weight: 232.23

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-Keto-D-gluconic acid potassium is a product generated by <i>Gluconobacter oxydans</i> . It can be effectively converted into L-(+)-tartaric acid.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.3061 mL	21.5304 mL	43.0608 mL
5 mM	0.8612 mL	4.3061 mL	8.6122 mL
10 mM	0.4306 mL	2.153 mL	4.3061 mL
50 mM	0.0861 mL	0.4306 mL	0.8612 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481