


Ethinyl estradiol 3-b-D-glucuronide

Chemical Properties

CAS No. :	60134-76-1
Formula:	C ₂₆ H ₃₂ O ₈
Molecular Weight:	472.53
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ethinyl estradiol 3-b-D-glucuronide is a type of biochemical reagent used for glycoscience research. Glycoscience focuses on the structure, synthesis, biology, and evolution of sugars, involving carbohydrate chemistry, glycans formation and degradation enzymology, protein-glycan recognition, and the role of glycans in biological systems. This field is closely related to basic research, biomedical sciences, and biotechnology.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1163 mL	10.5813 mL	21.1627 mL
5 mM	0.4233 mL	2.1163 mL	4.2325 mL
10 mM	0.2116 mL	1.0581 mL	2.1163 mL
50 mM	0.0423 mL	0.2116 mL	0.4233 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481