


Barium carbonate

Chemical Properties

CAS No. :	513-77-9
Formula:	CBaO3
Molecular Weight:	197.34
Storage:	Store at RT Actual storage temperature shall be subject to the COA.



Biological Description

Description	Barium carbonate,AR,99% is a reagent used in biochemical reactions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.0674 mL	25.337 mL	50.674 mL
5 mM	1.0135 mL	5.0674 mL	10.1348 mL
10 mM	0.5067 mL	2.5337 mL	5.0674 mL
50 mM	0.1013 mL	0.5067 mL	1.0135 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481