


Phenylmethyl β -D-arabinopyranoside

Chemical Properties

CAS No. :	5329-50-0
Formula:	C ₁₂ H ₁₆ O ₅
Molecular Weight:	240.25
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Phenylmethyl β -D-arabinopyranoside is a biochemical reagent utilized in glycobiology research. Glycobiology examines carbohydrate structure, synthesis, biology, and evolution. It encompasses carbohydrate chemistry, glycan formation and degradation enzymology, protein-glycan recognition, and the role of glycans in biological systems. This field is closely linked to fundamental research, biomedicine, and biotechnology.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.1623 mL	20.8117 mL	41.6233 mL
5 mM	0.8325 mL	4.1623 mL	8.3247 mL
10 mM	0.4162 mL	2.0812 mL	4.1623 mL
50 mM	0.0832 mL	0.4162 mL	0.8325 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481