


Phenylmethyl 2-(acetylamino)-2-deoxy-4,6-O-(phenylmethylene)- α -D-glucopyranoside

Chemical Properties

CAS No. :	13343-63-0
Formula:	C ₂₂ H ₂₅ N ₁ O ₆
Molecular Weight:	399.44
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Phenylmethyl 2-(acetylamino)-2-deoxy-4,6-O-(phenylmethylene)- α -D-glucopyranoside is a biochemical reagent used in glycobiology research. Glycobiology focuses on the study of carbohydrates, including their structure, synthesis, biology, and evolution. This field involves carbohydrate chemistry, glycan formation and degradation enzymology, protein-glycan recognition, and the role of glycans in biological systems, and is closely related to fundamental research, biomedicine, and biotechnology.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5035 mL	12.5175 mL	25.035 mL
5 mM	0.5007 mL	2.5035 mL	5.007 mL
10 mM	0.2504 mL	1.2518 mL	2.5035 mL
50 mM	0.0501 mL	0.2504 mL	0.5007 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481