


5,6-Isopropylidene-L-gulonic acid γ -lactone

Chemical Properties

CAS No. :	94697-68-4
Formula:	C ₉ H ₁₄ O ₆
Molecular Weight:	218.2
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	5,6-Isopropylidene-L-gulonic acid γ -lactone is a biochemical reagent used in the study of glycobiology. Glycobiology examines the structure, synthesis, biology, and evolution of sugars, encompassing carbohydrate chemistry, the enzymology of glycans' formation and degradation, protein-glycan interactions, and the role of glycans in biological systems. This field is closely linked to fundamental research, biomedicine, and biotechnology.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.583 mL	22.9148 mL	45.8295 mL
5 mM	0.9166 mL	4.583 mL	9.1659 mL
10 mM	0.4583 mL	2.2915 mL	4.583 mL
50 mM	0.0917 mL	0.4583 mL	0.9166 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481