


p-Nitrophenyl 2-acetamido-3,4,6-tri-O-acetyl-β-D-glucopyranoside

Chemical Properties

CAS No. : 13089-27-5
 Formula: C₂₀H₂₄N₂O₁₁
 Molecular Weight: 468.41
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	p-Nitrophenyl 2-acetamido-3,4,6-tri-O-acetyl-β-D-glucopyranoside has potential applications in glycobiology research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1349 mL	10.6744 mL	21.3488 mL
5 mM	0.427 mL	2.1349 mL	4.2698 mL
10 mM	0.2135 mL	1.0674 mL	2.1349 mL
50 mM	0.0427 mL	0.2135 mL	0.427 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481