


## Ethyl nonadecanoate

## Chemical Properties

CAS No. : 18281-04-4  
Formula: C<sub>21</sub>H<sub>42</sub>O<sub>2</sub>  
Molecular Weight: 326.56  
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Ethyl nonadecanoate is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0622 mL	15.3111 mL	30.6222 mL
5 mM	0.6124 mL	3.0622 mL	6.1244 mL
10 mM	0.3062 mL	1.5311 mL	3.0622 mL
50 mM	0.0612 mL	0.3062 mL	0.6124 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481