


2-Chloro-4-nitrophenyl-β-D-galactopyranoside

Chemical Properties

CAS No. : 123706-60-5
 Formula: C₁₂H₁₄ClNO₈
 Molecular Weight: 335.69
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Chloro-4-nitrophenyl-β-D-galactopyranoside (2-Chloro-4-nitrophenyl-beta-D-galactopyranoside) serves as a chromogenic substrate for assessing enzyme activity. It is commonly employed in molecular biology research to detect and measure β-galactosidase activity (biosynth: EC03318).
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9789 mL	14.8947 mL	29.7894 mL
5 mM	0.5958 mL	2.9789 mL	5.9579 mL
10 mM	0.2979 mL	1.4895 mL	2.9789 mL
50 mM	0.0596 mL	0.2979 mL	0.5958 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481