


(S)-1,2-Dilauroyl-sn-glycero-3-phosphocholine

Chemical Properties

CAS No. : 127641-86-5
 Formula: C32H64NO8P
 Molecular Weight: 621.83
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	(S)-1,2-Dilauroyl-sn-glycero-3-phosphocholine is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6082 mL	8.0408 mL	16.0816 mL
5 mM	0.3216 mL	1.6082 mL	3.2163 mL
10 mM	0.1608 mL	0.8041 mL	1.6082 mL
50 mM	0.0322 mL	0.1608 mL	0.3216 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481