


2-Nitrophenyl 2,3,4-Tri-O-acetyl- β -D-glucuronide,methyl ester

Chemical Properties

CAS No. : 55274-44-7
 Formula: C₁₉H₂₁NO₁₂
 Molecular Weight: 455.37
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Nitrophenyl 2,3,4-Tri-O-acetyl- β -D-glucuronide,methyl ester is a reagent used in biochemical reactions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.196 mL	10.9801 mL	21.9602 mL
5 mM	0.4392 mL	2.196 mL	4.392 mL
10 mM	0.2196 mL	1.098 mL	2.196 mL
50 mM	0.0439 mL	0.2196 mL	0.4392 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481