


Potassium D-tartrate monobasic

Chemical Properties

CAS No. : 57341-16-9
 Formula: C₄H₅KO₆
 Molecular Weight: 188.18
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Potassium D-tartrate monobasic is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.3141 mL	26.5703 mL	53.1406 mL
5 mM	1.0628 mL	5.3141 mL	10.6281 mL
10 mM	0.5314 mL	2.657 mL	5.3141 mL
50 mM	0.1063 mL	0.5314 mL	1.0628 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481