


## DL-Thyronine

## Chemical Properties

CAS No. :	1034-10-2
Formula:	C <sub>15</sub> H <sub>15</sub> NO <sub>4</sub>
Molecular Weight:	273.28
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	DL-Thyronine protects trypsin from photosensitized inactivation caused by riboflavin-5'-phosphate (FMN).
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6593 mL	18.2963 mL	36.5925 mL
5 mM	0.7319 mL	3.6593 mL	7.3185 mL
10 mM	0.3659 mL	1.8296 mL	3.6593 mL
50 mM	0.0732 mL	0.3659 mL	0.7319 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481