


## 1-Docosanoyl-sn-glycero-3-phosphocholine

## Chemical Properties

CAS No. :	125146-65-8
Formula:	C30H62NO7P
Molecular Weight:	579.79
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	1-Docosanoyl-sn-glycero-3-phosphocholine is a reagent used in biochemical reactions.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7248 mL	8.6238 mL	17.2476 mL
5 mM	0.345 mL	1.7248 mL	3.4495 mL
10 mM	0.1725 mL	0.8624 mL	1.7248 mL
50 mM	0.0345 mL	0.1725 mL	0.345 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481