


Glycochenodeoxycholic acid 3-sulfate disodium

Chemical Properties

CAS No. : 66874-09-7
 Formula: C₂₆H₄₁NNa₂O₈S
 Molecular Weight: 573.65
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Glycochenodeoxycholic acid 3-sulfate disodium is a reagent used in biochemical reactions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7432 mL	8.7161 mL	17.4322 mL
5 mM	0.3486 mL	1.7432 mL	3.4864 mL
10 mM	0.1743 mL	0.8716 mL	1.7432 mL
50 mM	0.0349 mL	0.1743 mL	0.3486 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481