


Phenylmethyl 3-O-(phenylmethyl)-4,6-O-(phenylmethylene)- α -D-mannopyranoside

Chemical Properties

CAS No. : 62774-16-7
 Formula: C₂₇H₂₈O₆
 Molecular Weight: 448.51
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Phenylmethyl 3-O-(phenylmethyl)-4,6-O-(phenylmethylene)- α -D-mannopyranoside is a type of biochemical reagent employed in glycobiology research. Glycobiology focuses on the structure, synthesis, biology, and evolution of sugars, encompassing carbohydrate chemistry, glycan formation and degradation enzymology, protein-glycan recognition, and the functions of glycans within biological systems. This field is closely linked to fundamental research, biomedicine, and biotechnology.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2296 mL	11.148 mL	22.296 mL
5 mM	0.4459 mL	2.2296 mL	4.4592 mL
10 mM	0.223 mL	1.1148 mL	2.2296 mL
50 mM	0.0446 mL	0.223 mL	0.4459 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481