


## p-Nitrophenyl 2-acetamido-2-deoxy-3-O-acetyl-β-D-glucopyranoside

## Chemical Properties

CAS No. : 23262-57-9  
 Formula: C<sub>16</sub>H<sub>20</sub>N<sub>2</sub>O<sub>9</sub>  
 Molecular Weight: 384.34  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	p-Nitrophenyl 2-acetamido-2-deoxy-3-O-acetyl-β-D-glucopyranoside is a reagent used in biochemical reactions.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6019 mL	13.0093 mL	26.0186 mL
5 mM	0.5204 mL	2.6019 mL	5.2037 mL
10 mM	0.2602 mL	1.3009 mL	2.6019 mL
50 mM	0.052 mL	0.2602 mL	0.5204 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481