

4-Fluorohippuric acid

Chemical Properties


CAS No. : 366-79-0

Formula: C₉H₈FN₃

Molecular Weight: 197.16

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	4-Fluorohippuric acid is a primary metabolite of fluoropyridine, showing potential in research for analgesic and anticonvulsant properties.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.072 mL	25.3601 mL	50.7202 mL
5 mM	1.0144 mL	5.072 mL	10.144 mL
10 mM	0.5072 mL	2.536 mL	5.072 mL
50 mM	0.1014 mL	0.5072 mL	1.0144 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481