


4,6-Dichloro-3-formylcoumarin

Chemical Properties

CAS No. :	51069-87-5
Formula:	C ₁₀ H ₄ Cl ₂ O ₃
Molecular Weight:	243.04
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	4,6-Dichloro-3-formylcoumarin is a reagent used in biochemical reactions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.1145 mL	20.5727 mL	41.1455 mL
5 mM	0.8229 mL	4.1145 mL	8.2291 mL
10 mM	0.4115 mL	2.0573 mL	4.1145 mL
50 mM	0.0823 mL	0.4115 mL	0.8229 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481