


Sulfosuccinic acid

Chemical Properties

CAS No. : 5138-18-1
 Formula: C₄H₆O₇S
 Molecular Weight: 198.15
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Sulfosuccinic acid is a reagent used in biochemical reactions.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.0467 mL	25.2334 mL	50.4668 mL
5 mM	1.0093 mL	5.0467 mL	10.0934 mL
10 mM	0.5047 mL	2.5233 mL	5.0467 mL
50 mM	0.1009 mL	0.5047 mL	1.0093 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481