


1-Heptanoyl-2-hydroxy-sn-glycero-3-phosphocholine

Chemical Properties

CAS No. : 160118-49-0
 Formula: C₁₅H₃₂N₁O₇P
 Molecular Weight: 369.39
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	1-Heptanoyl-2-hydroxy-sn-glycero-3-phosphocholine is a reagent used in biochemical reactions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7072 mL	13.5358 mL	27.0717 mL
5 mM	0.5414 mL	2.7072 mL	5.4143 mL
10 mM	0.2707 mL	1.3536 mL	2.7072 mL
50 mM	0.0541 mL	0.2707 mL	0.5414 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481