


1,3-Propanediyl bismethanethiosulfonate

Chemical Properties

CAS No. :	55-96-9
Formula:	C ₅ H ₁₂ O ₄ S ₄
Molecular Weight:	264.41
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	1,3-Propanediyl bismethanethiosulfonate is a reagent used in biochemical reactions.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.782 mL	18.910 mL	37.8201 mL
5 mM	0.7564 mL	3.782 mL	7.564 mL
10 mM	0.3782 mL	1.891 mL	3.782 mL
50 mM	0.0756 mL	0.3782 mL	0.7564 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481