


1,2-Ethanediy l bismethanethiosulfonate

Chemical Properties

CAS No. : 55-95-8
 Formula: C₄H₁₀O₄S₄
 Molecular Weight: 250.38
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	1,2-Ethanediy l bismethanethiosulfonate is a reagent used in biochemical reactions.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9939 mL	19.9696 mL	39.9393 mL
5 mM	0.7988 mL	3.9939 mL	7.9879 mL
10 mM	0.3994 mL	1.997 mL	3.9939 mL
50 mM	0.0799 mL	0.3994 mL	0.7988 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481