


## Sulfo-KMUS sodium salt

## Chemical Properties

CAS No. :	211236-68-9
Formula:	C19H25N2NaO9S
Molecular Weight:	480.47
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Sulfo-KMUS (sodium salt) can serve as a crosslinking agent.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0813 mL	10.4065 mL	20.813 mL
5 mM	0.4163 mL	2.0813 mL	4.1626 mL
10 mM	0.2081 mL	1.0406 mL	2.0813 mL
50 mM	0.0416 mL	0.2081 mL	0.4163 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481