


## (Rac)-1-Palmitoyl-glycero-3-phosphocholine

### Chemical Properties

CAS No. : 17364-18-0  
 Formula: C<sub>24</sub>H<sub>50</sub>NO<sub>7</sub>P  
 Molecular Weight: 495.63  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



### Biological Description

Description	(Rac)-1-Palmitoyl-glycero-3-phosphocholine is a reagent used in biochemical reactions.
Targets(IC50)	Others

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0176 mL	10.0882 mL	20.1763 mL
5 mM	0.4035 mL	2.0176 mL	4.0353 mL
10 mM	0.2018 mL	1.0088 mL	2.0176 mL
50 mM	0.0404 mL	0.2018 mL	0.4035 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481