


## 2,3-Dihydroxypropyl (8Z)-8-heptadecenoate

## Chemical Properties

CAS No. :	749266-30-6
Formula:	C <sub>20</sub> H <sub>38</sub> O <sub>4</sub>
Molecular Weight:	342.51
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2,3-Dihydroxypropyl (8Z)-8-heptadecenoate is a reagent used in biochemical reactions.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9196 mL	14.5981 mL	29.1962 mL
5 mM	0.5839 mL	2.9196 mL	5.8392 mL
10 mM	0.292 mL	1.4598 mL	2.9196 mL
50 mM	0.0584 mL	0.292 mL	0.5839 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481