


ABEI

Chemical Properties

CAS No. : 66612-29-1
 Formula: C₁₄H₂₀N₄O₂
 Molecular Weight: 276.33
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	ABEI is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6189 mL	18.0943 mL	36.1886 mL
5 mM	0.7238 mL	3.6189 mL	7.2377 mL
10 mM	0.3619 mL	1.8094 mL	3.6189 mL
50 mM	0.0724 mL	0.3619 mL	0.7238 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481