


4-Aminophenyl β-D-glucuronide

Chemical Properties

CAS No. : 21080-66-0
 Formula: C₁₂H₁₅NO₇
 Molecular Weight: 285.25
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	4-Aminophenyl β-D-glucuronide is a reagent used in biochemical reactions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5057 mL	17.5285 mL	35.057 mL
5 mM	0.7011 mL	3.5057 mL	7.0114 mL
10 mM	0.3506 mL	1.7528 mL	3.5057 mL
50 mM	0.0701 mL	0.3506 mL	0.7011 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481